
progetti[®]
Medical Equipment Solutions

RESCUE SAM

*Automated external defibrillator (AED)
Defibrillatore esterno automatico (DAE)*

DEFIBRILLATORS

Defibrillator

RESCUE SAM is a automated external defibrillator available in semi and fully-automated versions designed for providing simple and direct voice prompts and indications for a straightforward rescue operation.

It is lightweight and battery powered for maximum portability.

RESCUE SAM is designed to treat Ventricular Fibrillation (VF) and Fast Ventricular Tachycardia (FVT).

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
RESCUE SAM

TECHNICAL SPECIFICATIONS

Operating mode: semi-automated or fully-automated versions
Waveform: BTE (Biphasic Truncated Exponential); compensated by patient's impedance
Energy:
- adult: nominal 200 J with disposable pads for adults
- pediatric patients: nominal 50 J with optional disposable pediatric pads
Shock energy accuracy: +/-15%
Protocol: BLS-D with voice prompts and indicators guide protocol
Pacemaker: detection and rejection of artifacts
Charge control: automatic by software (charging control system once VF or FVT are detected), Analysis time: about 6 seconds
Charging time from "shock advised": <8 seconds, with new battery
Voice instructions: detailed voice messages indicate to the rescuer the actions needed
Controls: shock button, info-button, On/Off button
Indicators: LEDs (different colors), info-button
Charge complete indicator:
- The acoustic charging sound (beep) ceases
- Voice prompts: "press the red shock button" (semi-automated version) or "shock will be delivered in 3...2...1" (Fully-automated version)
- The shock button lights up steady red (semi-automated version only)
Disarm: once charged, the RESCUE SAM disarms if:
- The SHOCK button is not pressed within 15 seconds after being (for the semi-automated version only) or
- The Power ON/OFF button is pressed to turn OFF
Shock delivery: shock is delivered if the SHOCK button is pressed (semi-automated version) or automatically at the end of the countdown (fully-automated version)
Shock delivery configuration: via adult defibrillator pads in the anterior-anterior (Lead II) position or via reduced-energy pediatric pads in the anterior-posterior position
Applied parts type: BF-type defibrillation proof
Battery: 12V_{DC}, 4200mAh, LiMnO₂, non-rechargeable
Battery capacity: 200 shocks or 4 hours of continuous operating time
Self test:
- Power-On Self-Test: performed at every startup; it checks unit's basic functionality
- Scheduled Self-Tests: automatically run daily, weekly, and monthly; they verify hardware and software integrity
- Manual Self-Test: launched by removing and reinserting the battery pack; it tests the full system (charge and shock functions included)
Data recording and transmission:
- Free dedicated software to download data
- USB adapter for data transmission
- Data stored: up to 8 sessions; up to 60 minutes each. The device records event information, including ECG analysis and shockable/non-shockable rhythm classification
Ingress protection: IP 54
Size: 29 x 28 x 9 cm
Weight: 2.2 kg with battery pack

STANDARD SET:

- main unit
- disposable multifunction electrodes for adult patients
- non-rechargeable battery pack

OPTIONALS:

- disposable multifunction electrodes for adults
- disposable multifunction electrodes for pediatric patients
- transport bag
- wall mounted cabinet
- USB cable adapter for PC data management
- wall mount
- emergency kit

Discover Rescue Sam in action

- scan the QR code to watch the video presentation!

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PROGETTI S.r.l.